



**\*Price Available on Request**

**Hand Held Homogenizer**

A handheld homogenizer is a compact, portable laboratory device used to rapidly blend, mix, emulsify, suspend, or disrupt biological samples, plant/animal tissues, or other liquid/semi-solid substances. These devices use a high-speed motor to drive a rotor-stator mechanism, generating high shear forces to break down materials

Make	EASELAB
Capacity	10 - 250 ml
Standard Head	10 mm
Maximum Viscosity	5000 MPAS
Speed	5000 - 30000 rpm
Speed Mode	6 gears
Speed regulation	Step less speed regulation
Minimum Processing Capacity	10 ml
Noise	< 50 DB (A)
Protection Class	IP 30
Power	160 W
Power Supply	AC 220 v / 50 - 60 Hz
Packing Dimension (L x W x H)	480 x 280 x 110 mm
Gross Weight	3 kg

**For Order and other enquiry you may contact:**

<https://www.easelab.in>  
[easelab.india@gmail.com](mailto:easelab.india@gmail.com)  
 +91 730 984 4444