



***Price Available on Request**

Ice Flake Machine

An ice flake machine creates small, thin, super cooled ice pieces (flakes) by spraying water onto a chilled rotating drum or cylinder, which freezes the water into a thin sheet that a scraper then shaves off. These machines are widely used in food preservation (like seafood, produce, dairy), healthcare, and labs for cooling samples and products, providing ideal, fast-cooling ice for applications needing rapid temperature reduction and excellent contact.

Make	EASELAB
Model	ELIF-001
Capacity of Ice flake (kg per 24 hrs.)	30 Kg.
Storage Capacity (kg ~ app.)	12 Kg.
Storage bin	ABS Plastic
Cooling Method	Air
Ambient temperature	Below 25°C for optimum working
Construction	Stainless Steel
Noise level	≤52 dB
Water Supply	3/4 inch BSP
Refrigerant	R290A
Input Power	180W
Power supply	AC 220V,50/60Hz
Dimension LxWxH (mm)	610 x 470 x 780mm

For Order and other enquiry you may contact:

<https://www.easelab.in>

easelab.india@gmail.com

+91 730 984 4444