



***Price Available on Request**

Mini Dry Bath Incubator

EASELAB Mini Dry bath incubator features precise temperature control resolution in a compact design to give you reproducible results. This microprocessor controlled device offers you flexibility to accommodate a variety of interchangeable heating blocks for your versatile applications. Used in Immunoassays, Melting points, Enzyme reactions, Boiling points, Incubation, Activation of cultures, Laboratory procedures.

Make	EASELAB
Model	ELDBI-001
Temperature Range	RT+5 °C-100 °C
Temp. Setting Range	0 °C~100 °C
Temp. Stability@100°C	±0.5 °C
Temp. Display Accuracy	0.1 °C
Operation	Stirring & Heating
Module temperature stability	≤ ± 0.3 °C
Over Temperature Protection	Available
Heating Time	5 to 10min
Ambient temperature	10°C ~30°C
Cooling Speed	natural cooling
Power Supply	AC 220V / 50-60Hz

STANDARD ACCESSORIES : Main Instrument with Stand Block, ABS Protection cover, Power Cord

For Order and other enquiry you may contact:

<https://www.easelab.in>

easelab.india@gmail.com

+91 730 984 4444