



***Price Available on Request**



Multi Position Magnetic Stirrer

A multi-position magnetic stirrer is a laboratory instrument designed to mix multiple samples simultaneously on a single platform, enhancing efficiency and saving space. These units are available with various numbers of stirring points and may include heating capabilities.

Make	EASELAB
Plate Material	PET Anti-Corrosion Plate
Temperature Range	NA
Speed	50 - 1500 RPM
Control	Asynchronous Control (Synchronization can be set)
Number of Stirring Points	4
Max Stirring Volume	4 x 500 ml
Plate Size	100 x 110 mm (each)
Teflon Paddle Size	30 mm
Display	LCD
Time Setting	Yes
Power Supply	AC 220 W / 50-60 Hz
Ambient Temperature	5 ~ 40°C
Packing Dimension (LxWxH)	380 x 290 x 115mm
Gross Weight	2 kg

STANDARD ACCESSORIES

Main Instrument, Magnetic Stirring Bar, Power Cord, Instruction Manual, Test Certificate.

For Order and other enquiry you may contact:

<https://www.easelab.in>

easelab.india@gmail.com

+91 730 984 4444