



*Price Available on Request

Orbital Shaker

These shakers are designed to agitate samples in a circular motion, typically at high speeds. These units work best for bacterial cell growth, tissue dissociation, and other applications that don't require a lot of speed. These Shakers have a digital control panel for fine-tuning the mixing action.

Make	EASELAB
Model	ELOS-01
Display	LCD
Speed Range	40-200 RPM
Plate Type	Universal Tray Platform
Plate Size	250x250x50 mm
Orbital Motion	Circular
Max Load	2kg
Orbit diameter	20mm
Motor Type	Bush less Motor
Timing Range	Yes
Ambient Temperature	5°C to 40°C
Relative Humidity	80%
Power Supply	AC 220V/50-60Hz
Packing Dimension (LxWxH)	410 x 300 x 170mm
Gross Weight	5kg

Silent Feature:

- Universal Tray Platform for holding various containers.
- Suitable for Various Containers
- Start/Pause Switch for start or pause the operation.
- Timing knob can set the duration of operation.
- Digital Display

STANDARD ACCESSORIES: Main Instrument, Power cord, O Ring, Instruction Manual, Test Certificate

For Order and other enquiry you may contact:

<https://www.easelab.in>
easelab.india@gmail.com
+91 730 984 4444