



**\*Price Available on Request**

**Tube Roller Mixture**

A tube roller mixer is a laboratory instrument designed for the gentle and continuous mixing of samples in various tubes and cylindrical containers, widely used in clinical diagnostics, molecular biology, and biochemistry.

|                           |  |
|---------------------------|--|
| Make                      | EASELAB  |
| Speed                     | 20-80 rpm  |
| Ambient Temperature       | +5°C - +40°C   |
| Relative Humidity Range   | ≤80%   |
| Type                      | Rotating Shaking   |
| Max. Capacity             | 10 x 50 ml Centrifuge tube/0.5 kg  |
| Control Mode              | Digital key  |
| Display                   | LCD  |
| Time Setting              | YES  |
| Motor                     | Bushless   |
| Rollers                   | 6 Rollers  |
| Roller Size               | φ 6-34 mm  |
| Roller Length             | 280 mm   |
| Power Supply              | AC 220 V / 50 - 60 Hz  |
| Packing Dimension (LxWxH) | 510 x 305 x 145 mm   |
| Gross Weight              | 3 kg   |
| STANDARD ACCESSORIES      | Main Instrument, Power Cable, Power cord, Instruction Manual, Test Certificate |

**STANDARD ACCESSORIES:** Main Instrument, Power Cable , Power cord, Instruction Manual , Test Certificate

**For Order and other enquiry you may contact:**

<https://www.easelab.in>

[easelab.india@gmail.com](mailto:easelab.india@gmail.com)

+91 730 984 4444