



***Price Available on Request**

Ultrasonic Cleaner (Sonicator)

Ultrasonic Cleaners is a machinery which cleans surfaces by using high frequency sound waves. It can be used to clean all sorts of materials such as glass, plastics, metal, rubber as well as ceramics. A few of it's common applications are in hospitals, laboratories etc.

SALIENT FEATURES

- Normal and soft power control for the cleaning spots .
- Degas function for improve the cleaning effect & SS304
- Time Setting:1-99minutes adjustable.
- Heating Temperature RT+5oC ~ to 80oC
- User-friendly and clear panel
- SS304 Construction Inner Tank
- SS304 Construction Outer Tank
- SS304 Construction Lid (Top)
- Memory function (Test Setting)
- LID & Basket is Standard Supply

STANDARD ACCESSORIES : Main Instrument, Lid - SS Top, Power cable, Basket, Instruction Manual, Test Certificate

For Order and other enquiry you may contact:

<https://www.easelab.in>

easelab.india@gmail.com

+91 730 984 4444