



**\*Price Available on Request**

### **Vortex Mixture**

A **Vortex Mixer** (or Vortexer) is a common lab tool that rapidly mixes small vials of liquids by spinning a rubber cup, creating a swirling vortex in the liquid for tasks like reuse pending cells, mixing reagents for DNA extraction, or preparing blood samples, offering gentle to vigorous agitation for consistency in experiments.

Make	EASELAB
Model	ELVM-001
Speed Range	0 – 3000 rpm
Speed Timer	1~99 h 59 m 59 s
Display	White LCD
Shaking Movement	Orbital
Run Type	Continuous/Touch operation
Power	60W
Permissible ambient temperature	5-40°C
Permissible relative humidity	80%

#### **SILENT FEATURE:**

- Make EASELAB
- Operating Mode - 2 Types - Auto run sensor, continuous run used for various Mixing application
- Supply with Standard Head
- So silicone foot pads
- Shockproof, Bush less Motor
- Aluminum Alloy Base stable body

**For Order and other enquiry you may contact:**

<https://www.easelab.in>  
[easelab.india@gmail.com](mailto:easelab.india@gmail.com)  
 +91 730 984 4444