



***Price Available on Request**

PH Meter

Automatic Temperature Compensation. There is a microprocessor chip inside the meter and the appearance of meter looks nice with user friendly operation. Some of the features are as under, Large LCD display with blue background simultaneous display of pH and temperature, It has intelligent functions for calibration, data storage / Memory (upto 25 prints) Automatically display the percentage of slope after calibration is over. Equipped with electrode holder and calibration buffer solution, which is more convenient to use.

Make	EASELAB
Range	0 to 14.00 pH
pH Resolution	0.1 pH / 0.01 pH
Accuracy	± 0.01 pH
Calibration Points	up to 3 points
pH Buffer options	USA (pH 4.01, 7.00, 10.01)
mV Range	+1999 mV
Resolution & Accuracy	± 1 % FS + 1 digit & 1 mV
Temperature Range	0°C to 100°C
Accuracy	± 0.5°C
Resolution	0.1°C
Temperature Compensation	Automatic
Connector	BNC
Display	LCD
Memory	Upto 25 prints
Power Requirements	DC Adaptor, Using AC 220V,50/60Hz
Packing Dimensions (LxWxH)	300 × 250 × 130 mm
Weight	1kg

STANDARD ACCESSORIES : pH Meter, pH Electrode, pH Buffer - 4.01, 7.00, Power Adapter, Electrode Stand, ATC probe, Instruction Manual, Test Certificate.

For Order and other enquiry you may contact:

<https://www.easelab.in>

easelab.india@gmail.com

+91 730 984 4444