



**\*Price Available on Request**

### Dispensing Peristaltic Pump

A dispensing peristaltic pump is a high-precision, positive displacement pump that moves fluids by compressing flexible tubing with rotating rollers, ensuring the fluid touches only the tube, thus maintaining sterility. They are ideal for accurate, contamination-free dosing in pharmaceutical, lab, and food industries.

Make	EASELAB
Model	ELDPP-001
Speed Range	0.1~150rpm .
Speed resolution	0.1rpm (speed ≤100rpm )
	1rpm (speed >100rpm )
Timer	Yes
Communication interface	RS485 communication interface, MODBUS protocol is available
Standard Tube	18#
Power Consumption	<30W
Ambient Temperature	10°C ~ 40°C
Relative Humidity	<80%
IP grade	IP31
Circuit Board	Dust & Moisture proof
Power supply	AC 240V 50/60Hz
Housing	ABS Plastic
Pump Head / Channel	1
Brand	EASELAB
Application	Laboratory Use

**For Order and other enquiry you may contact:**

<https://www.easelab.in>  
[easelab.india@gmail.com](mailto:easelab.india@gmail.com)  
 +91 730 984 4444